

PATIENT EDUCATION

IGG Continuous Subcutaneous Infusion

Properly administering your medication at home is important to your safety. In the event of an emergency, call 911.

SUPPLIES:

- Freedom 60 Pump(s)
- Freedom 60/Precision Flow- tubing (**tubing to be changed DAILY**)
- Medication vials
- Sharps container
- (2) Sodium Chloride saline) prefilled syringes
- (1) Heparin prefilled syringe (**if needed, per your Plan of Treatment**)
- (1) Light Blue syringe tip cap, if indicated
- Alcohol/antiseptic wipes
- 50 cc Syringe(s) for medication
- Mini transfer device(s)
- Chloraprep swab sticks (3) **OR** Alcohol swabs if allergy or sensitivity to chloraprep
- Sub Q needle sets (25g to 27g) with transparent dressing included (single vs multiple needles per orders)

SASH FLUSH PROCEDURE:

- S Saline ___ ml
- A Administer medication as instructed by your nurse
- S Saline ___ ml
- H Heparin ___ ml (if required)

PROCEDURE:

1. Clean work area. Wash hands thoroughly for at least 20 seconds. Gather supplies and equipment. Remove medication from the refrigerator, as instructed by your nurse.
2. Check medication label for name, drug, frequency, and expiration. Inspect the medication for any cracks, leaks, particulate matter, and clarity of medication. Contact us for any discrepancies or concerns.
3. Draw up IGG per plan of treatment into the empty 50cc syringe, using a mini transfer device or needle.

PATIENT EDUCATION

IGG Continuous Subcutaneous Infusion

IMPORTANT SAFETY ALERT: Verify that you have the correct Freedom 60/Precision flow tubing. The rate printed on the package must match the rate ordered by your doctor (check prescription label).

4. Remove tubing from package.
5. Attach Freedom 60/Precision tubing to the luer disc end of medication syringe.
6. Place syringe into pump as instructed by your nurse, follow instruction sheet titled Medication Administration via Freedom 60 Pump.
7. Ensure the Freedom pump is in the **OFF** position and the black tab is at the end of its track. To move the black tab, wind the large knob clockwise until click is heard. **DO NOT OVERWIND.**

MEDICATION ADMINISTRATION:

8. Insert the medication syringe with tubing attached, plunger end first into the clear cover on the pump. Ensure luer disk end of tubing is firmly seated in the nose of the pump. (See instructional video)
9. Prime the Freedom 60/Precision Tubing by turning pump to **ON** position. Remove end cap UNTIL a drop of fluid is noted at end of tubing, reattach end cap. **Do NOT prime the needle. NOTE: Hold the syringe in place UNTIL the black tab engages with syringe to prevent disengagement.**
10. Turn the pump to **OFF** position.
11. Choose appropriate SQ site for patient. (Appropriate insertion sites include the abdomen, upper buttocks, lateral thigh, or lateral hip. Select a new site with each new needle placement.)
12. Cleanse site with a Chloraprep (or alcohol if chloraprep allergy) using a back-and-forth friction scrub for at least 30 seconds and let air dry for 60/Precision seconds. Allow the site to dry thoroughly. **Do not blot site dry.**
13. Carefully remove the protective cover from the SC needle, keeping the exposed needle sterile.
14. At the needle insertion site selected, pinch the skin, and insert the SQ needle at a 90-degree angle into the skin fold.

PATIENT EDUCATION

IGG Continuous Subcutaneous Infusion

15. Apply clear dressing over the needle insertion sites.
16. Secure excess tubing to skin using adhesive tape.
17. Remove protective cap from the end of Freedom 60/Precision Tubing and attach tubing to the injection cap on the end of your IV catheter. **Be careful to avoid touching the tip of the tubing.**
18. Turn pump ON to begin infusion. Pump will automatically begin to infuse medication. Administer medication via Freedom pump per your Plan of Treatment, monitoring infusion progress while medication is being dispensed. **Note:** Pump is silent during medication infusion, as instructed by your nurse.

AT THE END OF THE INFUSION:

19. When medication syringe is empty, Infusion is complete. Turn pump to **OFF** Position.
20. Wash hands thoroughly, for at least 20 seconds.
21. Gently remove transparent dressing while securing SQ needle.
22. Remove SQ needle(s) from the skin.
23. Rewind the knob clockwise, moving the black tab to the end of the track.
24. Remove the empty medication syringe from the pump. Leave the empty syringe attached to the tubing.
25. Discard the used syringe, tubing, and needle into the Sharps container.
26. Dispose of any IGG that is leftover in the vial, into the Sharps container.
27. Cover site with a gauze pad. You may apply direct pressure but **avoid rubbing site.**
28. Apply band aid(s).

PATIENT EDUCATION

IGG Continuous Subcutaneous Infusion

29. Wash hands thoroughly, for at least 20 seconds.

NOTE: Freedom 60/Precision tubing is to be disposed of daily.

****All supplies and waste can be double bagged and disposed of with your household garbage.****

NOTES:

CAREPATHrx

The information included in this document is intended for informational and educational purposes only. It is not intended to replace medical advice offered by a physician or other qualified health provider. If you think you are experiencing a medical emergency, call 911 immediately.

Copyright© 2021 by CarepathRx. All rights reserved.